


# CGRP ELISA kits



## Discover our CGRP ELISA kit offer for rat and human samples

- Ready to use validated kits (including quality control)
- High sensitivity with low background
- Detailed protocols with advice on sample preparation

# Calcitonin Gene Related Peptide (CGRP) biomarker: an exciting target for migraine & headache treatment. Obtain high specificity for human and rat samples.

- Sample volume :  
**100 µL** (25µL when extraction is performed)
- Limit of detection :  
**0.7 pg/mL** (in EIA buffer) and **2 pg/mL** (in plasma)
- Standard curve range :  
7.81 – 1000 pg/mL (**human**) / 3.91 – 500 pg/mL (**rat**)

## Specify for human kit

Human CGRP-I//II	100%
Human CGRP-a/β	100%
Rat CGRP-a/β	120%
Rat CGRP-I//II	120%

## Specify for rat kit

Rat CGRP-I//II	100%
Rat CGRP-a/β	100%
Human CGRP-a/β	83%
Human CGRP-I//II	83%

**CGRP (rat) ELISA kit : A05482**  
**CGRP (human) ELISA kit : A05481**

\* These kits are widely used across renowned Academic and Governmental Research Institutes, Biotechnology and Pharmaceutical companies

## VISIT OUR WEBSITE

Distributed By:



In Italy: **Vinci-Biochem Srl**

Via Ponte di Bagnolo, 10

50059 Vinci (Firenze) Italia

Tel: 0571 568 147 e 0571 568 135

[vb@vincibiochem.it](mailto:vb@vincibiochem.it) - <https://www.vincibiochem.it>



Bertin Technologies SAS

Parc d'activités du Pas du Lac

10 bis, avenue Ampère

78180 Montigny-le-Bretonneux, FRANCE

01 39 30 60 00

EU webstore: [bertin-bioreagent.com](http://bertin-bioreagent.com) | US webstore: [bertin-corp.com](http://bertin-corp.com)

[WWW.BERTIN-TECHNOLOGIES.COM](http://WWW.BERTIN-TECHNOLOGIES.COM)